

Smad1(Phospho-Ser206) Antibody

Catalog No: #11579

Package Size: #11579-1 50ul #11579-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Smad1(Phospho-Ser206) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of Smad1 only when phosphorylated at serine 206.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around phosphorylation site of serine 206 (P-H-S(p)-P-T) derived from Human Smad1.
Conjugates	Unconjugated
Target Name	Smad1
Modification	Phospho
Other Names	Mad-related protein 1; SMAD 1; hSMAD1; Transforming growth factor-beta-signaling protein 1;
Accession No.	Swiss-Prot: Q15797NCBI Protein: NP_001003688.1
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 60kd

Western blotting: 1:500~1:1000

Images



Western blot analysis of extracts from HeLa cells untreated or treated with PMA using Smad1(Phospho-Ser206) Antibody #11579.

Background

Transcriptional modulator activated by BMP (bone morphogenetic proteins) type 1 receptor kinase. SMAD1 is a receptor-regulated SMAD (R-SMAD).
Hiroaki Seto, et al. (2004) J Clin Invest. March 1; 113(5): 718

Note: This product is for in vitro research use only and is not intended for use in humans or animals.