

Acetyl-CoA Carboxylase (Phospho-Ser79) Antibody


Catalog No: #11584

Package Size: #11584-1 50ul #11584-2 100ul


Description

Product Name	Acetyl-CoA Carboxylase (Phospho-Ser79) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide-KLH
Conjugates	Unconjugated
Target Name	Acetyl-CoA Carboxylase
Modification	Phospho
Other Names	ACC, ACAC, ACC1, ACCA,
Accession No.	Swiss-Prot:P08581 Gene ID:4233
Concentration	1.0mg/ml

Images



Western blot analysis of extracts from 293 cells untreated or treated with AICAR using Acetyl-CoA Carboxylase (Phospho-Ser79) Antibody #11584.



Western blotting analysis using Acetyl-CoA Carboxylase (Phospho-Ser79) Antibody #11584.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.