

TNNI3 (Phospho-Ser43) Antibody

Catalog No: #12040

Package Size: #12040-1 50ul #12040-2 100ul

Description


Product Name	TNNI3 (Phospho-Ser43) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of TNNI3 only when phosphorylated at Serine 43.
Immunogen Type	Peptide-KLH
Immunogen Description	The antiserum was produced against synthesized peptide derived from human TNNI3 around the phosphorylation site of Ser43.
Conjugates	Unconjugated
Target Name	TNNI3
Modification	Phospho
Other Names	CMH7, RCM1, cTnI, CMD2A, TNNC1
Accession No.	Swiss-Prot#: P19429; NCBI Gene#: 7137; NCBI Protein#: NP_000354.4
SDS-PAGE MW	26kd
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Predicted MW: 26kd

Western blotting: 1:500~1:1000

Images



Western blot analysis of lysates from Jurkat cells, using TNNI3 (Phospho-Ser43) Antibody. The lane on the right is blocked with the phospho peptide.



Immunohistochemistry analysis of paraffin-embedded human heart, using TNNI3 (Phospho-Ser43) Antibody. The picture on the right is blocked with the phospho peptide.

Background

Troponin I (TnI), along with troponin T (TnT) and troponin C (TnC), is one of 3 subunits that form the troponin complex of the thin filaments of striated muscle. TnI is the inhibitory subunit; blocking actin-myosin interactions and thereby mediating striated muscle relaxation. The TnI subfamily contains three genes: TnI-skeletal-fast-twitch, TnI-skeletal-slow-twitch, and TnI-cardiac. This gene encodes the TnI-cardiac protein and is exclusively expressed in cardiac muscle tissues. Mutations in this gene cause familial hypertrophic cardiomyopathy type 7 (CMH7) and familial restrictive cardiomyopathy (RCM). [provided by RefSeq, Jul 2008].

Note: This product is for in vitro research use only and is not intended for use in humans or animals.