

Pan Phospho-Serine/Threonine Rabbit pAb

Catalog No: #12098

Package Size: #12098-1 50ul #12098-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Pan Phospho-Serine/Threonine Rabbit pAb
Host Species	rabbit
Clonality	polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse;Rat
Immunogen Type	Synthetic peptide.
Conjugates	Unconjugated
SDS-PAGE MW	>10kDa
Formulation	PBS with 0.09% Sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:1000


Images



Western blot analysis of lysates from NIH/3T3 cells using Pan Phospho-Serine/Threonine Rabbit pAb at 1:400 dilution.



Western blot analysis of lysates from HeLa cells using Pan Phospho-Serine/Threonine Rabbit pAb at 1:400 dilution.



Western blot analysis of lysates from C6 cells using Pan Phospho-Serine/Threonine Rabbit pAb at 1:400 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.