

HNRPA0 (Phospho-Ser84) Antibody

Catalog No: #12686

Package Size: #12686-1 50ul #12686-2 100ul


Description

Product Name	HNRPA0 (Phospho-Ser84) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu
Specificity	HNRPA0 (Phospho-Ser84) Antibody detects endogenous levels of HNRPA0 only when phosphorylated at Ser84
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human HNRPA0 (Phospho-Ser84)
Conjugates	Unconjugated
Target Name	HNRPA0
Modification	Phospho
Other Names	HNRNPA0, HnRNA binding protein, HnRNP A0, HNRPA0
Accession No.	Swiss-Prot#: Q13151NCBI Gene ID: 10949
Target Species	human
Calculated MW	31kd
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western Blot analysis of lysates of 1. HeLa cell and 2. HeLa cell treated with 100ng/mL LPS for 30min, using primary antibody at 1:1000 dilution.

Published Papers

el at., A tumor-specific modulation of heterogeneous ribonucleoprotein A0 promotes excessive mitosis and growth in colorectal cancer cells.In Cell Death Dis on 2020 Apr 17 by Konishi H, Fujiya M,et al..PMID:32303675, , (2020)

[PMID:32303675](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.