

LATS2 (Phospho-Ser872) Antibody

Catalog No: #12875

Package Size: #12875-1 50ul #12875-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LATS2 (Phospho-Ser872) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Phospho-LATS2 (S872) Antibody detects endogenous levels of LATS2 only when phosphorylated at S872
Immunogen Type	Peptide-KLH
Immunogen Description	A synthesized peptide derived from human LATS2 (Phospho-Ser872)
Conjugates	Unconjugated
Other Names	FLJ13161 antibody Kinase phosphorylated during mitosis protein antibody KPM antibody Large tumor suppressor homolog 2 antibody Large tumor suppressor homolog 2 Drosophila antibody Large tumor suppressor kinase 2 antibody LATS large tumor suppressor Drosophila homolog 2 antibody LATS large tumor suppressor homolog 2 antibody Lats2 antibody LATS2_HUMAN antibody Serine threonine kinase kpm antibody Serine threonine protein kinase kpm antibody Serine threonine-protein kinase kpm antibody Serine threonine-protein kinase LATS2 antibody Warts like kinase antibody Warts-like kinase antibody
Accession No.	Swiss-Prot#:Q9NRM7 NCBI Gene ID26524
Calculated MW	150
Concentration	1.0mg mL
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+) pH 7.4 150mM NaCl 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB dilution:1:1000

Published Papers

et al., Identification of the kinase STK25 as an upstream activator of LATS signalling. In Nat Commun on 2019 Apr 4 by Lim S, Hermance N, et al..PMID:30948712, , (2019)

[PMID:30948712](#)

et al., The CalcR-PKA-Yap1 Axis Is Critical for Maintaining Quiescence In Muscle Stem Cells. In Cell Rep on 2019 Nov 19 by Zhang L, Noguchi YT, et al..PMID:31747590, , (2019)

[PMID:31747590](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.