

STING (Phospho-Ser366) Antibody

Catalog No: #13260

Package Size: #13260 100ul


Description

Product Name	STING (Phospho-Ser366) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human; Mouse; Rat
Specificity	This antibody detects endogenous levels of Human/Rat STING (phospho-Ser366) or mouse STING (phospho-Ser365)
Immunogen Description	Synthesized phosho peptide around human STING (Ser366)
Conjugates	Unconjugated
Target Name	STING
Other Names	endoplasmic reticulum IFN stimulator; Endoplasmic reticulum interferon stimulator; ERIS; FLJ38577; hMITA; hSTING; Mediator of IRF3 activation; MITA; Mitochondrial mediator of IRF3 activation; MPYS; N terminal methionine proline tyrosine serine plasma membrane tetraspanner; NET23; Stimulator of interferon genes; Stimulator of interferon genes protein; STING; TM173_HUMAN; Tmem173; Transmembrane protein 173;
Accession No.	Swiss-Prot#:Q86WV6NCBI Gene ID:340061
Calculated MW	38kDa
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20 °C.Stable for 12 months from date of receipt

Application Details

WB 1:1000-1:2000

Images



Published Papers

Byers Lauren A;Cardnell Robert J;Cargill Kasey R;Gay Carl M;Kundu Kiran;Ramkumar Kavya;Shen Li;Stewart C Allison;Wang Jing;Zhang Bingnan et al., SLFN11 biomarker status predicts response to lurbinectedin as a single agent and in combination with ATR inhibition in small cell lung cancer, , (2021)

[PMID:35004241](#)

Bishu Wang;Xingguo Tang;Chuntao Xiao;Zhijie Yu;Huaben Bo;Jie Wang;Jinquan Wang et al., Nucleus-targeted ruthenium(II) complex triggers immunogenic cell death and sensitizes melanoma to anti-PD-1 therapy by activating cGAS-STING pathway., , (2025)

[PMID:40022761](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.