

PP2A(Phospho-Y307) Rabbit mAb

Catalog No: #13369

Package Size: #13369-1 50ul #13369-2 100ul

Description

Product Name	PP2A(Phospho-Y307) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	ST49-05
Purification	ProA affinity purified
Applications	WB, ICC/IF, IHC
Species Reactivity	Hu, Ms, Rt
Immunogen Description	Synthetic phospho-peptide corresponding to residues surrounding Tyr307 of human PP2A.
Conjugates	Unconjugated
Other Names	PP2A A antibody PP2A alpha antibody PP2A B antibody PP2A beta antibody PP2A-alpha antibody PP2AA_HUMAN antibody PP2Ac antibody PP2CB antibody PPP2CA antibody PPP2CB antibody Protein phosphatase 2 catalytic subunit alpha isoform antibody Protein phosphatase 2 catalytic subunit beta isoform antibody Replication protein C antibody RP C antibody RP-C antibody Serine/threonine protein phosphatase 2A catalytic subunit alpha isoform antibody Serine/threonine protein phosphatase 2A catalytic subunit beta isoform antibody Serine/threonine-protein phosphatase 2A catalytic subunit alpha isoform antibody
Accession No.	Swiss-Prot#:P62714
Calculated MW	Predicted band size: 35 kDa
SDS-PAGE MW	Observed band size: 35 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C


Application Details

WB: 1:500-1:2000

ICC/IF: 1:50-1:200

IHC: 1:50-1:200

Images



All lanes : PP2A(Phospho-Y307) Rabbit mAb at 1/1k dilution

Lane 1 : A431 whole cell lysates

Lane 2 : F9 whole cell lysates

Lane 3 : PC12 whole cell lysates

Lane 4 : Mouse spleen lysates

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution

Predicted band size: 35 kDa

Observed band size: 35 kDa

Exposure time: 12 seconds



All lanes : PP2A alpha (Phospho-Tyr307) Rabbit mAb at 1/1k dilution

Lane 1 : Rat Kidney lysates

Lane 2 : Mouse Liver lysates

Lane 3 : A431 treated with 100ng/ml EGF for 20min whole cell lysates

Lane 4 : A431 whole cell lysates

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution


Predicted band size: 35 kDa

Observed band size: 35 kDa


Exposure time: 10 seconds



Formalin-fixed, paraffin-embedded human colon tissue stained for PP2A(Phospho-Y307) using 13369 at 1/100 dilution in immunohistochemical analysis.



Formalin-fixed, paraffin-embedded human lung carcinoma tissue stained for PP2A(Phospho-Y307) using 13369 at 1/100 dilution in immunohistochemical analysis.



Immunocytochemistry/Immunofluorescence
 PP2A(Phospho-Y307) antibody (13369)
 ICC/IF staining of PP2A(Phospho-Y307) in Hela cells. Cells were fixed with 4% Paraformaldehyde permeabilized with 0.1% Triton X-100. Samples were incubated with 13369(red) at a working dilution of 1/100.HSP70 Monoclonal antibody (green) was diluted at 1:200. Goat Anti Rabbit Alexa Fluor 647 was diluted at 1:1000. Goat Anti Mouse Alexa Fluor 488 was diluted at 1:1000.

Background

Protein phosphatase type 2A (PP2A) is an essential protein serine/threonine phosphatase that is conserved in all eukaryotes. PP2A is a key enzyme within various signal transduction pathways as it regulates fundamental cellular activities such as DNA replication, transcription, translation, metabolism, cell cycle progression, cell division, apoptosis and development. The core enzyme consists of catalytic C and regulatory A (or PR65) subunits, with each subunit represented by α and β isoforms. Additional regulatory subunits belong to four different families of unrelated proteins. Both the B (or PR55) and B' regulatory protein families contain α , β , γ and δ isoforms, with the B' family also including an ϵ protein. B'' family proteins include PR72, PR130, PR59 and PR48 isoforms, while striatin (PR110) and SG2NA (PR93) are both members of the B''' regulatory protein family. These B subunits competitively bind to a shared binding site on the core A subunit. This variable array of holoenzyme components, particularly regulatory B subunits, allows PP2A to act in a diverse set of functions. PP2A function is regulated by expression, localization, holoenzyme composition and post-translational modification. Phosphorylation of PP2A at Tyr307 by Src occurs in response to EGF or insulin and results in a substantial reduction of PP2A activity. Reversible methylation on the carboxyl group of Leu309 of PP2A has been observed. Methylation alters the conformation of PP2A, as well as its localization and association with B regulatory subunits.

References

Published Papers

Lei Chen;Xia Zhao;Rui Sheng;Philip Lazarovici;Wenhua Zheng et al., Artemisinin alleviates astrocyte overactivation and neuroinflammation by modulating the IRE1/NF- κ B signaling pathway in in vitro and in vivo Alzheimer's disease models., (2025)

[PMID:39826816](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.