

TCTP (phospho-Ser46) rabbit pAb

Catalog No: #13487

Package Size: #13487-1 50ul #13487-2 100ul

Description

Product Name	TCTP (phospho-Ser46) rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	This antibody detects endogenous levels of Human Mouse Rat TCTP (phospho-Ser46)
Immunogen Description	Synthesized phosho peptide around human TCTP (Ser46)
Conjugates	Unconjugated
Other Names	Translationally-controlled tumor protein (TCTP) (Fortilin) (Histamine-releasing factor) (HRF) (p23)
Accession No.	Swiss Prot:P13693GenelD:7178
SDS-PAGE MW	24
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.