

RSK2 (phospho-Ser227) rabbit pAb

Catalog No: #13542

Package Size: #13542-1 50ul #13542-2 100ul

Description

Product Name	RSK2 (phospho-Ser227) rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	This antibody detects endogenous levels of Human Mouse Rat RSK2 (phospho-Ser227)
Immunogen Description	Synthesized phospho peptide around human RSK2 (Ser227)
Conjugates	Unconjugated
Accession No.	Swiss Prot:Q15418
SDS-PAGE MW	85
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.