

## RIP (phospho-Ser166) rabbit pAb

Catalog No: #13557

Package Size: #13557-1 50ul #13557-2 100ul

## Description

|                       |                                                                                                                                                                                 |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | RIP (phospho-Ser166) rabbit pAb                                                                                                                                                 |
| Host Species          | Rabbit                                                                                                                                                                          |
| Clonality             | Polyclonal                                                                                                                                                                      |
| Purification          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                                       |
| Applications          | WB                                                                                                                                                                              |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                 |
| Specificity           | This antibody detects endogenous levels of Human RIP (phospho-Ser166)                                                                                                           |
| Immunogen Description | Synthesized phosho peptide around human RIP (Ser166)                                                                                                                            |
| Conjugates            | Unconjugated                                                                                                                                                                    |
| Other Names           | Receptor-interacting serine/threonine-protein kinase 1 (EC 2.7.11.1) (Cell death protein RIP)<br>(Receptor-interacting protein 1) (RIP-1) (Serine/threonine-protein kinase RIP) |
| Accession No.         | Swiss Prot:Q13546GenelD:8737                                                                                                                                                    |
| SDS-PAGE MW           | 73                                                                                                                                                                              |
| Concentration         | 1 mg/ml                                                                                                                                                                         |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                         |
| Storage               | -20°C/1                                                                                                                                                                         |

## Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.