

# Kv4.2 (phospho Ser616) Polyclonal Antibody

Catalog No: #13756

Package Size: #13756-1 50ul #13756-2 100ul

## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | Kv4.2 (phospho Ser616) Polyclonal Antibody                                                                            |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | IHC;IF;ELISA                                                                                                          |
| Species Reactivity    | Human,Mouse,Rat                                                                                                       |
| Specificity           | Phospho-Kv4.2 (S616) Polyclonal Antibody detects endogenous levels of Kv4.2 protein only when phosphorylated at S616. |
| Immunogen Description | Synthesized phospho-peptide around the phosphorylation site of human Kv4.2 (phospho Ser616)                           |
| Conjugates            | Unconjugated                                                                                                          |
| Other Names           | KCND2; KIAA1044; Potassium voltage-gated channel subfamily D member 2; Voltage-gated potassium channel subunit Kv4.2  |
| Accession No.         | Swiss Prot:Q9NZV8GeneID:3751                                                                                          |
| Concentration         | 1 mg/ml                                                                                                               |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | -20°C/1                                                                                                               |

## Application Details

IHC 1:100-1:300; ELISA 1:5000; IF 1:50-200

## Background

potassium voltage-gated channel subfamily D member 2(KCND2) Homo sapiens Voltage-gated potassium (Kv) channels represent the most complex class of voltage-gated ion channels from both functional and structural standpoints. Their diverse functions include regulating neurotransmitter release, heart rate, insulin secretion, neuronal excitability, epithelial electrolyte transport, smooth muscle contraction, and cell volume. Four sequence-related potassium channel genes - shaker, shaw, shab, and shal - have been identified in Drosophila, and each has been shown to have human homolog(s). This gene encodes a member of the potassium channel, voltage-gated, shal-related subfamily, members of which form voltage-activated A-type potassium ion channels and are prominent in the repolarization phase of the action potential. This member mediates a rapidly inactivating, A-type outward potassium current which is not under the control of the N terminus as i

Note: This product is for in vitro research use only and is not intended for use in humans or animals.