

JNK1/2/3 (phospho Thr183/Y185) Polyclonal Antibody

Catalog No: #13769

Package Size: #13769-1 50ul #13769-2 100ul

Description

Product Name	JNK1/2/3 (phospho Thr183/Y185) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IF/ICC, WB, IHC-p, ELISA
Species Reactivity	Human, Mouse, Rat, Chicken (tested by our customer)
Specificity	Phospho-JNK1/2/3 (T183/Y185) Polyclonal Antibody detects endogenous levels of JNK1/2/3 protein only when phosphorylated at T183/Y185.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human JNK1/2/3 around the phosphorylation site of Thr183 and Tyr185. AA range: 151-200
Conjugates	Unconjugated
Other Names	MAPK8; JNK1; PRKM8; SAPK1; SAPK1C; Mitogen-activated protein kinase 8; MAP kinase 8; MAPK 8; JNK-46; Stress-activated protein kinase 1c; SAPK1c; Stress-activated protein kinase JNK1; c-Jun N-terminal kinase 1; MAPK9; JNK2; PRKM9; SAPK1A; Mi
Accession No.	Swiss Prot: P45983/P45984/P53779 GeneID: 5599/5601/5602
SDS-PAGE MW	46+54
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

IF: 1:50-200 WB 1:500-2000, IHC 1:50-300 IHC 1:50-300

Background

mitogen-activated protein kinase 8 (MAPK8) Homo sapiens The protein encoded by this gene is a member of the MAP kinase family. MAP kinases act as an integration point for multiple biochemical signals, and are involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation and development. This kinase is activated by various cell stimuli, and targets specific transcription factors, and thus mediates immediate-early gene expression in response to cell stimuli. The activation of this kinase by tumor-necrosis factor alpha (TNF-alpha) is found to be required for TNF-alpha induced apoptosis. This kinase is also involved in UV radiation induced apoptosis, which is thought to be related to cytochrome c-mediated cell death pathway. Studies of the mouse counterpart of this gene suggested that this kinase play a key role in T cell proliferation, apoptosis and differentiation. Several alternatively spliced

Note: This product is for in vitro research use only and is not intended for use in humans or animals.