

# Catenin $\delta$ -1 (phospho-Tyr904) rabbit pAb

Catalog No: #14006

Package Size: #14006-1 50ul #14006-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Catenin $\delta$ -1 (phospho-Tyr904) rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human Catenin $\delta$ -1 (phospho-Tyr904)
Immunogen Description	Synthesized phosho peptide around human Catenin $\delta$ -1 (Tyr904)
Conjugates	Unconjugated
Other Names	Catenin delta-1 (Cadherin-associated Src substrate) (CAS) (p120 catenin) (p120(ctn)) (p120(cas))
Accession No.	Swiss Prot:O60716GeneID:1500
SDS-PAGE MW	108
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

## Application Details

WB 1:1000-2000

## Background

catenin delta 1 (CTNND1) Homo sapiens This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and alternative splicing result in many different isoforms being translated. Not all of the full-length nature of the described transcript variants have been determined. Read-through transcription also exists between this gene and the neighboring upstream thioredoxin-related transmembrane protein 2 (TMX2) gene. [provided by RefSeq, Dec 2010],

Note: This product is for in vitro research use only and is not intended for use in humans or animals.