

## Atg14 (phospho-Ser29) rabbit pAb

Catalog No: #14061

Package Size: #14061-1 50ul #14061-2 100ul

## Description

|                       |                                                                                                                   |
|-----------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name          | Atg14 (phospho-Ser29) rabbit pAb                                                                                  |
| Host Species          | Rabbit                                                                                                            |
| Clonality             | Polyclonal                                                                                                        |
| Purification          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.         |
| Applications          | WB                                                                                                                |
| Species Reactivity    | Human, Mouse                                                                                                      |
| Specificity           | This antibody detects endogenous levels of Human Mouse Atg14 (phospho-Ser29)                                      |
| Immunogen Description | Synthesized phospho peptide around human Atg14 (Ser29)                                                            |
| Conjugates            | Unconjugated                                                                                                      |
| Other Names           | Beclin 1-associated autophagy-related key regulator (Barkor) (Autophagy-related protein 14-like protein) (Atg14L) |
| Accession No.         | Swiss Prot:Q6ZNE5GeneID:22863                                                                                     |
| SDS-PAGE MW           | 55                                                                                                                |
| Concentration         | 1 mg/ml                                                                                                           |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                           |
| Storage               | -20°C/1                                                                                                           |

## Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.