

## 53BP1 (phospho-Ser1618) rabbit pAb

Catalog No: #14115

Package Size: #14115-1 50ul #14115-2 100ul

## Description

|                       |                                                                                                           |
|-----------------------|-----------------------------------------------------------------------------------------------------------|
| Product Name          | 53BP1 (phospho-Ser1618) rabbit pAb                                                                        |
| Host Species          | Rabbit                                                                                                    |
| Clonality             | Polyclonal                                                                                                |
| Purification          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Applications          | WB                                                                                                        |
| Species Reactivity    | Human,Rat                                                                                                 |
| Specificity           | This antibody detects endogenous levels of Human Rat 53BP1 (phospho-Ser1618)                              |
| Immunogen Description | Synthesized phosho peptide around human 53BP1 (Ser1618)                                                   |
| Conjugates            | Unconjugated                                                                                              |
| Other Names           | Tumor suppressor p53-binding protein 1 (53BP1) (p53-binding protein 1) (p53BP1)                           |
| Accession No.         | Swiss Prot:Q12888GeneID:7158                                                                              |
| SDS-PAGE MW           | 213                                                                                                       |
| Concentration         | 1 mg/ml                                                                                                   |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                   |
| Storage               | -20°C/1                                                                                                   |

## Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.