

# PP2A alpha (Phospho-Tyr307) Rabbit mAb

Catalog No: #14129

Package Size: #14129-1 50ul #14129-2 100ul


## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | PP2A alpha (Phospho-Tyr307) Rabbit mAb                                                             |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB IHC ICC/IF IP                                                                                   |
| Species Reactivity    | Human Mouse Rat                                                                                    |
| Specificity           | Phospho-PP2A alpha (Y307) Antibody detects endogenous levels of total Phospho-PP2A alpha (Y307)    |
| Immunogen Description | A synthesized peptide derived from human Phospho-PP2A alpha (Y307)                                 |
| Conjugates            | Unconjugated                                                                                       |
| Other Names           | eplication protein C; PP2A-alpha; Replication protein C; RP-C; PPP2CA; MGC786                      |
| Accession No.         | Uniprot:P67775                                                                                     |
| Calculated MW         | 35kDa                                                                                              |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.                   |

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50

## Images



Western blot analysis of Phospho-PP2A alpha (Y307) in (1) Rat kidney lysate; (2) A431 cell lysate treated with EGF.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.