

ATM (Phospho-Ser1981) Rabbit mAb

Catalog No: #14136

Package Size: #14136-1 50ul #14136-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ATM (Phospho-Ser1981) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human
Specificity	Phospho-ATM (S1981) Antibody detects endogenous levels of total Phospho-ATM (S1981)
Immunogen Description	A synthesized peptide derived from human Phospho-ATM (S1981)
Conjugates	Unconjugated
Other Names	kinase ATM; Serine-protein kinase ATM
Accession No.	Uniprot:Q13315
Calculated MW	370kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of Phospho-ATM (Ser1981) in (1) HEK293 cell lysate; (2) HEK293 cell lysate treated with Doxorubicin.



Immunohistochemical analysis of paraffin-embedded human liver, using Phospho-ATM (S1981) Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.