

DNA PKcs (Phospho-Ser2056) Rabbit mAb

Catalog No: #14139

Package Size: #14139-1 50ul #14139-2 100ul

Description

Product Name	DNA PKcs (Phospho-Ser2056) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	Phospho-DNA PKcs (S2056) Antibody detects endogenous levels of total Phospho-DNA PKcs (S2056)
Immunogen Description	A synthesized peptide derived from human Phospho-DNA PKcs (S2056)
Conjugates	Unconjugated
Other Names	DNA- PKcs, DNA-dependent protein kinase catalytic subunit, DNPK1, EC 2.7.11.1, P460, PRKD, PRKDC, XRCC7, kinase DNA-PK
Accession No.	Uniprot:P78527
Calculated MW	469kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

IHC:1:50~1:200

ICC/IF:1:50~1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.