

p53 (Phospho-Ser392) Rabbit mAb

Catalog No: #14142

Package Size: #14142-1 50ul #14142-2 100ul


Description

Product Name	p53 (Phospho-Ser392) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	Phospho-p53 (S392) Antibody detects endogenous levels of total Phospho-p53 (S392)
Immunogen Description	A synthesized peptide derived from human Phospho-p53 (S392)
Conjugates	Unconjugated
Other Names	Antigen NY-CO-13, Cellular tumor antigen p53, Phosphoprotein p53, TP53, Tumor suppressor p53
Accession No.	Uniprot:P04637
Calculated MW	53kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of Phospho-p53 (Ser392) expression in HEK293 whole cell lysate.



Immunohistochemical analysis of paraffin-embedded human breast, using Phospho-p53 (S392) Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.