

EGFR (Phospho-Tyr1068) Rabbit mAb

Catalog No: #14144

Package Size: #14144-1 50ul #14144-2 100ul


Description

Product Name	EGFR (Phospho-Tyr1068) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse
Specificity	Phospho-EGFR (Y1068) Antibody detects endogenous levels of total Phospho-EGFR (Y1068)
Immunogen Description	A synthesized peptide derived from human Phospho-EGFR (Y1068)
Conjugates	Unconjugated
Other Names	ERBB1, Epidermal growth factor receptor precursor; Receptor protein-tyrosine kinase ErbB-1; kinase EGFR; HER1; EGFR; SA7; NISBD2;
Accession No.	Uniprot:P00533
Calculated MW	170kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of EGFR phosphorylation expression in A431 cell treated with EGF.



Immunohistochemical analysis of paraffin-embedded human breast cancer, using Phospho-EGFR (Y1068) Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.