

CDK1/2/3 (Phospho-Thr14) Rabbit mAb

Catalog No: #14153

Package Size: #14153-1 50ul #14153-2 100ul


Description

Product Name	CDK1/2/3 (Phospho-Thr14) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	Phospho-CDK1/2/3 (T14) Antibody detects endogenous levels of total Phospho-CDK1/2/3 (T14)
Immunogen Description	A synthesized peptide derived from human Phospho-CDK1/2/3 (T14)
Conjugates	Unconjugated
Other Names	CDK1; CDK2; CDK3;
Accession No.	Uniprot:P06493/P24941/Q00526
Calculated MW	34kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of Phospho-CDK1/2/3 (T14) expression in (1) A431 cell lysate; (2) HeLa cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.