

Src (Phospho-Tyr419) Rabbit mAb

Catalog No: #14159

Package Size: #14159-1 50ul #14159-2 100ul

Description

Product Name	Src (Phospho-Tyr419) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	Phospho-Src (Y419) Antibody detects endogenous levels of total Phospho-Src (Y419)
Immunogen Description	A synthesized peptide derived from human Phospho-Src (Y419)
Conjugates	Unconjugated
Other Names	p60-Src; c-Src; pp60c-src; Avian sarcoma virus; Proto-oncogene c-Src; SRC Oncogene;
Accession No.	Uniprot:P12931
Calculated MW	60kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000


FC:

30

Images



Western blot analysis of Phospho-Src (Y419) expression in A431 cell lysate treated with pervanadate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.