

p53 (Phospho-Thr55) Rabbit mAb

Catalog No: #14163

Package Size: #14163-1 50ul #14163-2 100ul


Description

Product Name	p53 (Phospho-Thr55) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Phospho-p53 (T55) Antibody detects endogenous levels of total Phospho-p53 (T55)
Immunogen Description	A synthesized peptide derived from human Phospho-p53 (T55)
Conjugates	Unconjugated
Other Names	Antigen NY-CO-13; TP53; Tumor suppressor p53; Cellular tumor antigen p53; Phosphoprotein p53
Accession No.	Uniprot:P04637
Calculated MW	53kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of Phospho-p53 (T55) in HEK293 cell lysates treated with Calyculin A and Okadaic Acid.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.