

Retinoblastoma (Phospho-Ser807) Rabbit mAb

Catalog No: #14165

Package Size: #14165-1 50ul #14165-2 100ul


Description

Product Name	Retinoblastoma (Phospho-Ser807) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Phospho-Retinoblastoma (S807) Antibody detects endogenous levels of total Phospho-Retinoblastoma (S807)
Immunogen Description	A synthesized peptide derived from human Phospho-Retinoblastoma (S807)
Conjugates	Unconjugated
Other Names	P105-RB; PP105; PP110; RB-1; RB1; Retinoblastoma-associated protein;
Accession No.	Uniprot:P06400
Calculated MW	106kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of Phospho-Retinoblastoma (S807) expression in K562 cell lysate.



Immunohistochemical analysis of paraffin-embedded mouse spleen, using Phospho-Retinoblastoma (S807) Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.