

c-Myc (Phospho-Ser62) Rabbit mAb

Catalog No: #14168

Package Size: #14168-1 50ul #14168-2 100ul


Description

Product Name	c-Myc (Phospho-Ser62) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	Phospho-c-Myc (S62) Antibody detects endogenous levels of total Phospho-c-Myc (S62)
Immunogen Description	A synthesized peptide derived from human Phospho-c-Myc (S62)
Conjugates	Unconjugated
Other Names	Myc proto-oncogene protein; c-myc; c Myc; Myc;
Accession No.	Uniprot:P01106
Calculated MW	57kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:100~1:500IP:1:50

Images



Western blot analysis of Phospho-c-Myc (S62) expression in HeLa cell lysate.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma, using Phospho-c-Myc (S62) Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.