

## BTK (Phospho-Tyr223) Rabbit mAb

Catalog No: #14170

Package Size: #14170-1 50ul #14170-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                                                                                |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | BTK (Phospho-Tyr223) Rabbit mAb                                                                                                                |
| Host Species          | Rabbit                                                                                                                                         |
| Clonality             | Monoclonal                                                                                                                                     |
| Isotype               | Rabbit IgG                                                                                                                                     |
| Purification          | Affinity-chromatography                                                                                                                        |
| Applications          | WB IP                                                                                                                                          |
| Species Reactivity    | Human                                                                                                                                          |
| Specificity           | Phospho-BTK (Y223) Antibody detects endogenous levels of total Phospho-BTK (Y223)                                                              |
| Immunogen Description | A synthesized peptide derived from human Phospho-BTK (Y223)                                                                                    |
| Conjugates            | Unconjugated                                                                                                                                   |
| Other Names           | Agammaglobulinaemia tyrosine kinase; AGMX1; ATK; B cell progenitor kinase; BPK; Bruton's tyrosine kinase; EC 2.7.10.2; kinase Btk; Kinase EMB; |
| Accession No.         | Uniprot:Q06187                                                                                                                                 |
| Calculated MW         | 76kDa                                                                                                                                          |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.                                             |
| Storage               | Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.                                             |

## Application Details

WB:1:1000~1:2000IP:1:20

## Images



Western blot analysis of Phospho-BTK (Y223) expression in Raji cell lysate treated with pervanadate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.