

Stat5 (Phospho-Tyr694) Rabbit mAb

Catalog No: #14172

Package Size: #14172-1 50ul #14172-2 100ul


Description

Product Name	Stat5 (Phospho-Tyr694) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse
Specificity	Phospho-Stat5 (Y694) Antibody detects endogenous levels of total Phospho-Stat5 (Y694)
Immunogen Description	A synthesized peptide derived from human Phospho-Stat5 (Y694)
Conjugates	Unconjugated
Other Names	Signal transducer and activator of transcription 5A; STAT5; STAT5A; STAT5B; Transcription factor STAT5A; Transcription factor STAT5B; MGF;
Accession No.	Uniprot:P42229/P51692
Calculated MW	90kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of Phospho-Stat5 (Y694) expression in A431 cell lysate treated with EGF.



Immunohistochemical analysis of paraffin-embedded human thyroid, using Phospho-Stat5 (Y694) Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.