

EGFR (Phospho-Tyr1110) Rabbit mAb

Catalog No: #14176

Package Size: #14176-1 50ul #14176-2 100ul


Description

Product Name	EGFR (Phospho-Tyr1110) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR1196
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Phospho-EGFR (Y1086) Antibody detects endogenous levels of Phospho-EGFR (Y1086)
Immunogen Description	A synthesized peptide derived from human Phospho-EGFR (Y1086)
Conjugates	Unconjugated
Other Names	EC 2.7.10.1; ERBB1; Epidermal growth factor receptor precursor; Receptor protein-tyrosine kinase ErbB-1; kinase EGFR;
Accession No.	Uniprot:P00533
Calculated MW	Predicted band size: 134 kDa
SDS-PAGE MW	Observed band size: 150 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $^{\circ}$ C short term. Store at -20 $^{\circ}$ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500-1:2000

Images



All lanes : EGFR (Phospho-Tyr1110) Rabbit mAb at 1/1k dilution
Lane 1 : A431 cell lysate
Lane 2 : A431 100ng/ml EGF 10min cell lysate
lysate/proteins at 20 μ g per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 134 kDa
Observed band size: 150 kDa
Exposure time: 7 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.