

TrkB (Phospho-Tyr817) Rabbit mAb

Catalog No: #14177

Package Size: #14177-1 50ul #14177-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TrkB (Phospho-Tyr817) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	Phospho-TrkB (Y817) Antibody detects endogenous levels of Phospho-TrkB (Y817)
Immunogen Description	A synthesized peptide derived from human Phospho-TrkB (Y817)
Conjugates	Unconjugated
Other Names	NTRK2; gp145-TrkB; Trk-B; TRKB; Tropomyosin-related kinase B; TrkB tyrosine kinase; Tyrosine kinase receptor B;
Accession No.	Uniprot:Q16620
Calculated MW	140kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:1000IHC:1:50~1:100ICC/IF:1:50~1:100IP:1:30

Images



Western blot analysis of Phospho-TrkB (Y817) expression in SH-SY5Y cell lysate treated with BDNF.



Immunohistochemical analysis of paraffin-embedded mouse brain cancer, using Phospho-TrkB (Y817) Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.