

AMPK alpha 1/2 (Phospho-Thr183/Thr172) Rabbit mAb

Catalog No: #14179

Package Size: #14179-1 50ul #14179-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	AMPK alpha 1/2 (Phospho-Thr183/Thr172) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR2077
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB,IHC
Species Reactivity	Human Mouse Rat
Specificity	Phospho-AMPK alpha (T172) Antibody detects endogenous levels of Phospho-AMPK alpha (T172)
Immunogen Description	A synthesized peptide derived from human AMPK alpha
Conjugates	Unconjugated
Other Names	ACACA kinase; Acetyl-CoA carboxylase kinase; AMPK alpha 2 chain; AMPK subunit alpha-2; AMPK2; AMPKalpha2; HMGCR kinase; PRKAA; PRKAA2;
Accession No.	Uniprot:Q13131/P54646
Calculated MW	Predicted band size: 64,62 kDa
SDS-PAGE MW	Observed band size: 62 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200

Images



All lanes : AMPK alpha 1/2 (Phospho-Thr183/Thr172) Rabbit mAb at 1/1k dilution
 Lane 1 : Mouse spleen lysate
 Lane 2 : 293 cell lysate
 Lane 3 : A735 cell lysate
 Lane 4 : Mouse heart lysate
 lysate/proteins at 20 μ g per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 64,62 kDa
 Observed band size: 62 kDa
 Exposure time: 14 seconds



Formalin-fixed, paraffin-embedded human brain tissue stained for AMPK alpha (Phospho-Thr172) using 14179 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.