

ATF2 (Phospho-Thr71) Rabbit mAb

Catalog No: #14183

Package Size: #14183-1 50ul #14183-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ATF2 (Phospho-Thr71) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human
Specificity	Phospho-ATF2 (T71) Antibody detects endogenous levels of Phospho-ATF2 (T71)
Immunogen Description	A synthesized peptide derived from human Phospho-ATF2 (T71)
Conjugates	Unconjugated
Other Names	ATF2, Activating 2, CREB2, CREBP1, Cyclic-AMP-dependent ATF-2, HB16, MXBP protein, cAMP response element binding protein CRE-BP1;
Accession No.	Uniprot:P15336
Calculated MW	70kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200IP:1:30

Images



Western blot analysis of Phospho-ATF2 (T71) expression in (1) HeLa cell lysate; (2) HeLa cell lysate treated with Anisomycin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.