

SHP2 (Phospho-Tyr542) Rabbit mAb

Catalog No: #14185

Package Size: #14185-1 50ul #14185-2 100ul


Description

Product Name	SHP2 (Phospho-Tyr542) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human Mouse
Specificity	Phospho-SHP2 (Y542) Antibody detects endogenous levels of Phospho-SHP2 (Y542)
Immunogen Description	A synthesized peptide derived from human Phospho-SHP2 (Y542)
Conjugates	Unconjugated
Other Names	BPTP3; CFC; MGC14433; NS1; PTN11; PTP-1D; PTP-2C; PTP2C; PTPN11; SH-PTP2; SH-PTP3; SHP-2; Shp2; SHPTP2;
Accession No.	Uniprot:Q06124
Calculated MW	68kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of Phospho-SHP2 (Y542) expression in (1) NIH/3T3 cell lysates; (2) NIH/3T3 cell lysates treated with PDGF.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.