

Synapsin I (Phospho-Ser9) Rabbit mAb

Catalog No: #14187

Package Size: #14187-1 50ul #14187-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	Synapsin I (Phospho-Ser9) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR2212
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB,IHC
Species Reactivity	Human Mouse Rat
Specificity	Phospho-Synapsin I (S9) Antibody detects endogenous levels of Phospho-Synapsin I (S9)
Immunogen Description	A synthesized peptide derived from human Phospho-Synapsin I (S9)
Conjugates	Unconjugated
Other Names	Brain protein 4.1; SYN-1; synapsin I;
Accession No.	Uniprot:P17600
Calculated MW	Predicted band size: 74 kDa
SDS-PAGE MW	Observed band size: 77 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200

Images



All lanes : Synapsin I (Phospho-Ser9) Rabbit mAb at 1/1k dilution
Lane 1 : Human brain lysate
Lane 2 : Human brain lysate treated with AP lysate/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 74 kDa Observed band size: 77 kDa
Exposure time: 5 seconds



Formalin-fixed, paraffin-embedded human brain tissue stained for Synapsin I (Phospho-Ser9) using 14187 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.