

MSK1 (Phospho-Ser376) Rabbit mAb

Catalog No: #14188

Package Size: #14188-1 50ul #14188-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MSK1 (Phospho-Ser376) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Rat
Specificity	Phospho-MSK1 (S376) Antibody detects endogenous levels of Phospho-MSK1 (S376)
Immunogen Description	A synthesized peptide derived from human Phospho-MSK1 (S376)
Conjugates	Unconjugated
Other Names	90 kDa ribosomal protein S6 kinase polypeptide 5; EC 2.7.11.1; KS6A5; MSPK1; Nuclear mitogen- and stress-activated protein kinase-1; RLPK; RPS6KA5; kinase MSK1; ribosomal protein S6 kinase;
Accession No.	Uniprot:O75582
Calculated MW	90kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF1:50~1:200IP:1:50

Images



Western blot analysis of Phospho-MSK1 (S376) expression in HEK293 cell lysate treated with EGF.



Immunohistochemical analysis of paraffin-embedded human lung cancer, using Phospho-MSK1 (S376) Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.