

RPA2 (Phospho-Thr21) Rabbit mAb

Catalog No: #14189

Package Size: #14189-1 50ul #14189-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RPA2 (Phospho-Thr21) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	Phospho-RPA2 (T21) Antibody detects endogenous levels of Phospho-RPA2 (T21)
Immunogen Description	A synthesized peptide derived from human Phospho-RPA2 (T21)
Conjugates	Unconjugated
Other Names	REPA2; RF-A; RP-A; RPA2; RPA32; replication factor-A protein 2; replication protein A 32 kDa subunit; replication protein A2;
Accession No.	Uniprot:P15927
Calculated MW	32kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:5000~1:10000IP:1:20

Images



Western blot analysis of Phospho-RPA2 (T21) expression in HeLa cell lysate treated with Calyculin A.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.