

SMC1 (Phospho-Ser957) Rabbit mAb

Catalog No: #14190

Package Size: #14190-1 50ul #14190-2 100ul


Description

Product Name	SMC1 (Phospho-Ser957) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	Phospho-SMC1 (S957) Antibody detects endogenous levels of Phospho-SMC1 (S957)
Immunogen Description	A synthesized peptide derived from human Phospho-SMC1 (S957)
Conjugates	Unconjugated
Other Names	SB1.8/DXS423E protein; SM1A; SMC1A; SMC1L1; SMC1alpha protein; Sb1.8;
Accession No.	Uniprot:Q14683
Calculated MW	160kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of Phospho-SMC1 (S957) expression in (1) HeLa cell lysate; (2) HeLa cells lysate treated with AP.



Immunohistochemical analysis of paraffin-embedded human bladder cancer, using Phospho-SMC1 (S957) Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.