

HSL (Phospho-Ser853) Rabbit mAb

Catalog No: #14191

Package Size: #14191-1 50ul #14191-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HSL (Phospho-Ser853) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Phospho-HSL (S853) Antibody detects endogenous levels of Phospho-HSL (S853)
Immunogen Description	A synthesized peptide derived from human Phospho-HSL (S853)
Conjugates	Unconjugated
Other Names	Hormone sensitive lipase; LIPE; LIPS; HSL;
Accession No.	Uniprot:Q05469
Calculated MW	83kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:5000~1:10000

Images



Western blot analysis of Phospho-HSL (S853) expression in (1) Mouse muscle lysate; (2) Mouse muscle lysate treated with AP.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.