

Chk1 (Phospho-Ser296) Rabbit mAb

Catalog No: #14195

Package Size: #14195-1 50ul #14195-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Chk1 (Phospho-Ser296) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Phospho-Chk1 (S296) Antibody detects endogenous levels of Phospho-Chk1 (S296)
Immunogen Description	A synthesized peptide derived from human Phospho-Chk1 (S296)
Conjugates	Unconjugated
Other Names	Checkpoint, S. pombe, homolog of, 1; CHEK1; CHK1; CHK1 checkpoint homolog (S. pombe); Serine/threonine-protein kinase Chk1;
Accession No.	Uniprot:O14757
Calculated MW	54kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of Phospho-Chk1 (S296) expression in HEK293 cell lysate Treated with Calyculin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.