

YAP1 (Phospho-Ser127) Rabbit mAb

Catalog No: #14199

Package Size: #14199-1 50ul #14199-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	YAP1 (Phospho-Ser127) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	Phospho-YAP1 (S127) Antibody detects endogenous levels of Phospho-YAP1 (S127)
Immunogen Description	A synthesized peptide derived from human YAP1
Conjugates	Unconjugated
Other Names	65 kDa Yes associated protein; YAP; YAP 1; YAP2; YAP 65; YKI; Yorkie homolog;
Accession No.	Uniprot:P46937
Calculated MW	65kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:5000~1:20000IHC:1:50~1:200

Images



Western blot analysis of Phospho-YAP1 (S127) expression in (1) HeLa cell lysate; (2) HeLa cell lysate treated with FBS+calyculin A.



Immunohistochemical analysis of paraffin-embedded human stomach, using Phospho-YAP1 (S127) Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.