

EGFR (Phospho-Tyr1068) Rabbit mAb

Catalog No: #14205

Package Size: #14205-1 50ul #14205-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EGFR (Phospho-Tyr1068) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human
Specificity	Phospho-EGFR (Y1068) Antibody detects endogenous levels of Phospho-EGFR (Y1068)
Immunogen Description	A synthesized peptide derived from human EGFR
Conjugates	Unconjugated
Other Names	Epidermal growth factor receptor; Proto-oncogene c-ErbB-1; Receptor tyrosine-protein kinase erbB-1; ERBB; ERBB1; HER1; EGFR; Urogastrone;
Accession No.	Uniprot:P00533
Calculated MW	175kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:2000|ICC/IF:1:50~1:200|IP:1:100|FC:1:40

Images



Western blot analysis of Phospho-EGFR (Y1068) expression in (1) A431 cell lysate; (2) A431 cell treated with EGF lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.