

## PLK1 (Phospho-Thr210) Rabbit mAb

Catalog No: #14206

Package Size: #14206-1 50ul #14206-2 100ul


## Description

Product Name	PLK1 (Phospho-Thr210) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Phospho-PLK1 (T210) Antibody detects endogenous levels of Phospho-PLK1 (T210)
Immunogen Description	A synthesized peptide derived from human PLK1
Conjugates	Unconjugated
Other Names	Cell cycle regulated protein kinase; PLK 1; PLK; Polo like kinase 1 ; Serine/threonine protein kinase 13; Serine/threonine protein kinase PLK1; STPK 13;
Accession No.	Uniprot:P53350
Calculated MW	68kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:1000~1:2000

## Images



Western blot analysis of Phospho-PLK1 (T210) expression in (1) HeLa cell lysate; (2) HeLa cell treated with nocodazole lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.