

Paxillin (Phospho-Tyr118) Rabbit mAb

Catalog No: #14217

Package Size: #14217-1 50ul #14217-2 100ul


Description

Product Name	Paxillin (Phospho-Tyr118) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR2124
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Phospho-Paxillin (Y118) Antibody detects endogenous levels of Phospho-Paxillin (Y118)
Immunogen Description	A synthesized peptide derived from human Paxillin
Conjugates	Unconjugated
Other Names	Paired box protein Pax 1; PAX 1; PAX1; Paxillin alpha; Paxillin; PXN; PXN protein;
Accession No.	Uniprot:P49023
Calculated MW	Predicted band size: 65 kDa
SDS-PAGE MW	Observed band size: 68 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4/°C short term. Store at -20/°C long term. Avoid freeze / thaw cycle.


Application Details

WB: 1:500-1:2000

Images



All lanes : Paxillin (Phospho-Tyr118) Rabbit mAb at 1/1k dilution
Lane 1 : HeLa cell lysate
Lane 2 : HeLa cell treated with pervanadate. lysate/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 65 kDa Observed band size: 68 kDa
Exposure time: 6 seconds



All lanes : Paxillin (Phospho-Tyr118) Rabbit mAb at 1/1k dilution
Lane 1 : Hela
Lane 2 : Hela KD
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 65 kDa
Observed band size: 68 kDa
Exposure time: 4 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.