

AMPK alpha 2 (Phospho-Ser345) Rabbit mAb

Catalog No: #14218

Package Size: #14218-1 50ul #14218-2 100ul


Description

Product Name	AMPK alpha 2 (Phospho-Ser345) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Phospho-AMPK alpha 2 (S345) Antibody detects endogenous levels of Phospho-AMPK alpha 2 (S345)
Immunogen Description	A synthesized peptide derived from human AMPK alpha 2
Conjugates	Unconjugated
Other Names	ACACA kinase; Acetyl-CoA carboxylase kinase; AMPK subunit alpha-2; AMPK2; AMPKalpha2; HMGCR kinase; Hydroxymethylglutaryl-CoA reductase kinase; PRKAA; PRKAA2;
Accession No.	Uniprot:P54646
Calculated MW	64kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:2000

Images



Western blot analysis of Phospho-AMPK alpha 2 (S345) expression in (1) 293T cell lysate treated with Lambda Phosphatase; (2) 293T cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.