

SHIP (Phospho-Tyr1022) Rabbit mAb

Catalog No: #14220

Package Size: #14220-1 50ul #14220-2 100ul


Description

Product Name	SHIP (Phospho-Tyr1022) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human
Specificity	Phospho-SHIP (Tyr1022) Antibody detects endogenous levels of Phospho-SHIP (Tyr1022)
Immunogen Description	A synthesized peptide derived from human SHIP
Conjugates	Unconjugated
Other Names	Phosphatidylinositol 3,4,5-trisphosphate 5-phosphatase 1; Inositol polyphosphate-5-phosphatase of 145 kDa; SIP-145; SH2 domain-containing inositol 5'-phosphatase 1; SHIP-1; p150Ship; hp51CN; INPP5D; SHIP;
Accession No.	Uniprot:Q92835
Calculated MW	145kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:2000ICC/IF:1:50~1:200IP:1:20

Images



Western blot analysis of Phospho-SHIP (Y1020) expression in (1) Raji cell lysate; (2) Raji cell treated with pervanadate lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.