

AKT1 (Phospho-Ser129) Rabbit mAb

Catalog No: #14222

Package Size: #14222-1 50ul #14222-2 100ul

Description

Product Name	AKT1 (Phospho-Ser129) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Phospho-AKT1 (S129) Antibody detects endogenous levels of Phospho-AKT1 (S129)
Immunogen Description	A synthesized peptide derived from human AKT1
Conjugates	Unconjugated
Other Names	AKT 1; PKB; PKB-ALPHA; PRKBA; Protein Kinase B Alpha; Protein kinase B; Proto-oncogene c-Akt; RAC Alpha; RAC-alpha serine/threonine-protein kinase; RAC-PK-alpha;
Accession No.	Uniprot:P31749
Calculated MW	60kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.