

Cdc6 (Phospho-Ser54) Rabbit mAb

Catalog No: #14225

Package Size: #14225-1 50ul #14225-2 100ul


Description

Product Name	Cdc6 (Phospho-Ser54) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human
Specificity	Phospho-Cdc6 (S54) Antibody detects endogenous levels of total Phospho-Cdc6 (S54)
Immunogen Description	A synthesized peptide derived from human Phospho-Cdc6 (S54)
Conjugates	Unconjugated
Other Names	CDC6; CDC18L; CDC6-related protein; Cdc18-related protein; HsCDC6; HsCDC18; p62(cdc6);
Accession No.	Uniprot:Q99741
Calculated MW	63kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:1000ICC/IF:1:50~1:200

Images



Western blot analysis of Phospho-Cdc6 (S54) expression in (1) Raji cell lysate; (2) Raji + FBS cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.