

Raf1 (Phospho-Ser621) Rabbit mAb

Catalog No: #14227

Package Size: #14227-1 50ul #14227-2 100ul

Description

Product Name	Raf1 (Phospho-Ser621) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Phospho-Raf1 (S621) Antibody detects endogenous levels of total Phospho-Raf1 (S621)
Immunogen Description	A synthesized peptide derived from human Phospho-Raf1 (S621)
Conjugates	Unconjugated
Other Names	C-Raf; CRAF; Protein kinase raf 1; RAF1; Raf proto-oncogene serine/threonine-protein kinase; kinase Raf1;
Accession No.	Uniprot:P04049
Calculated MW	73kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.