

Rac1/Cdc42 (Phospho-Ser71) Rabbit mAb

Catalog No: #14229

Package Size: #14229-1 50ul #14229-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Rac1/Cdc42 (Phospho-Ser71) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Phospho-Rac1/Cdc42 (Ser71) Antibody detects endogenous levels of total Phospho-Rac1/Cdc42 (Ser71)
Immunogen Description	A synthesized peptide derived from human Phospho-Rac1/Cdc42 (Ser71)
Conjugates	Unconjugated
Other Names	CDC42; CDC42Hs; G25K; TKS; MIG5; Ras like protein TC25;
Accession No.	Uniprot:P60953/P63000
Calculated MW	24kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:1000

Images



Western blot analysis of Phospho-Rac1/Cdc42 (Ser71) expression in (1) A431 cell lysate treated with LP; (2) A431 cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.