

Hsp27 (Phospho-Ser78) Rabbit mAb

Catalog No: #14230

Package Size: #14230-1 50ul #14230-2 100ul


Description

Product Name	Hsp27 (Phospho-Ser78) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human
Specificity	Phospho-Hsp27 (S78) Antibody detects endogenous levels of total Phospho-Hsp27 (S78)
Immunogen Description	A synthesized peptide derived from human Phospho-Hsp27 (S78)
Conjugates	Unconjugated
Other Names	28 kDa heat shock protein; Estrogen-regulated 24 kDa protein; Growth-related 25 kDa protein; Heat shock 27 kDa protein; HSP 27; HSP25; HSPB1; P25; SRP27; Stress-responsive protein 27;
Accession No.	Uniprot:P04792
Calculated MW	27kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of Phospho-Hsp27 (S78) expression in (1) A431 cell lysate; (2) A431 cell lysate treated with Anisomycin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.